

Maria Cifuentes Girard ¹
Patrick Knight ²
Roger Giles ²
Gérard Hopfgartner ¹

1. University of Geneva
Life Sciences Mass Spectrometry
Geneva, Switzerland

2. Shimadzu Research Laboratory
Manchester, United Kingdom

References

1. Shvartsburg, A. A., A. Haris, R. Andrzejewski, A. Entwistle and R. Giles (2018). "Differential Ion Mobility Separations in the Low-Pressure Regime." *Analytical Chemistry* **90**(1): 936-943.

2. Method description BfR-PA-Tee-2.0/2014; Determination of pyrrolizidine alkaloids (PA) in plant material using SPE-LC-MS/MS

LCMS-vDMS8060 is intended for Research Use Only and is not for use in diagnostic procedures. This presentation may contain references to products that are not available in your country. All rights reserved. Information subject to change without notice.

Overview

- Separation of diastereomeric pyrrolizidine alkaloids using differential ion mobility spectrometry-mass spectrometry (DMS-MS) at low pressure (33 mbar).
- Combination of vacuum DMS-MS with a short column in trap/elute setup to increase throughput in quantitative analysis.
- Comparison of analytical performance of the short LC-vDMS-MS (5 min) versus a reported LC-MS assay (12 min) [2].

Introduction

In recent years, the monitoring of phytotoxic compounds such as, pyrrolizidine alkaloids (PAs) in food matrices has gained interest. PAs are a complex mixture of diastereomers and more than 600 and their N-oxides forms were identified in over 6000 plants. Many LC-MS and GC-MS assays have been developed, but these methods require extensive sample preparation and long analyses times. Differential ion mobility is powerful for the separation of isomers and for filtering chemical noise. In the present work, the potential of vacuum DMS (vDMS) for enhanced selectivity and short LC columns are investigated for the analysis of diastereomeric PAs in tea samples for improved sample throughput.

Four set of PA diastereomers were investigated. Intermedine, echinatine, lycopsamine, indicine (*m/z* 300) and their N-oxide forms; intermedine-N-oxide, echinatine-N-oxide, indicine-N-oxide, lycopsamine-N-oxide (*m/z* 316) have no selective fragments under CID and coeluted in LC. Senecivernine and senecionine (*m/z* 336) cannot be distinguished by MS/MS but they partially coeluted in LC. Jacobine is the only diastereomer of *m/z* 352 that has a selective fragment (*m/z* 352>155) and can be distinguished from retrorsine, senecivernine-N-oxide and senecionine-N-oxide (*m/z* 352); which have different retention times under isocratic conditions at 15% MeOH. The orthogonality of vDMS as second dimension to LC-MS allowed the separation of co-eluting analytes without compromising analysis time. Several parameters in vDMS were optimized to achieved separation of diastereomers such as nature of solvent, amplitude of the waveform, compensation voltage, temperature and pressure. Sample preparation was simplified due to the combination of short column and trap/elute setup to vDMS-MS detection. Finally, PA diastereomers were analyzed in tea samples using the new developed short LC-DMS-MS method and compared to a LC-MS reference method.

Analytes

PAs	LC-MS MS (DMS mode)					Short LC-DMS-MS (MRM mode)				
	Q1 (m/z)	Q3 (m/z)	CE (V)	IS	RT (min)	RT (min)	CV (V)	DV (V)		
1 Intermedine (InT)	300.1	94.1	27	Alp	4.79	0.93	-18.16	800		
2 Echinatine (Ech)	300.1	94.1	20	Alp	4.86	0.91	-19.46	800		
3 Lycopsamine (Lyc)	300.1	94.1	28	Alp	4.89	0.93	-18.06	800		
4 Indicine (InD)	300.1	94.1	25	Alp	4.87	0.92	-16.85	800		
5 Intermedine-N-oxide (InT-N)	316.1	111.1	30	Alp	5.34	1.30	-18.36	800		
6 Echinatine-N-oxide (Ech-N)	316.1	111.1	41	Alp	5.25	1.31	-22.47	800		
7 Lycopsamine-N-oxide (Lyc-N)	316.1	111.1	25	Alp	5.47	1.39	-21.07	800		
8 Indicine-N-oxide (InD-N)	316.1	111.1	30	Alp	5.36	1.30	-18.36	800		
9 Jacobine (Jac)	352.1	120.1	20	Epi	4.76	0.93	-13.34	800		
10 Retrorsine (Ret)	352.1	120.1	30	Epi	5.88	1.49	-12.54	800		
11 Senecivernine-N-oxide (SeV-N)	352.1	94.1	30	Epi	7.15	2.49	-14.45	800		
12 Senecionine-N-oxide (SeC-N)	352.1	94.1	30	Epi	7.28	2.82	-14.85	800		
13 Senecivernine (SeV)	336.1	120.1	25	Epi	6.88	2.23	-13.24	800		
14 Senecionine (SeC)	336.1	120.1	25	Epi	6.97	2.39	-13.34	800		
15 Atropine (Atp)	290.3	93.1	35	N.A.	6.53	1.95	-16.22	800		
16 Epi-jacobine (Epi)	352.1	155.1	20	N.A.	6.35	1.71	-14.15	800		




Figure 1. Schematic representation of the vacuum Differential Ion Mobility-Mass Spectrometry (vDMS-MS).[1]

Materials and methods

Sample preparation

Figure 2. Tea sample preparation workflow: A) follow the BFR protocol [2] that includes extraction with 0.05 M sulfuric acid by sonication, centrifugation, neutralization, filtration, solid phase extraction, evaporation and reconstitution. B) a simplified protocol that excluded the solid phase extraction, evaporation and reconstitution steps for LC-DMS-MS analysis in trap/elute mode.

Short LC-DMS-MS method





Figure 3. Short LC-DMS-MS configuration. A trap column Reprosil –Pur C18-AQ (10x 0.5 mm, 5 µm), and a short column Kinetex XB-C18 (50 x 1.0 mm, 5 µm) was used to analyzed the samples in 2 steps: front-flush injection (a) and back-flush elution of analytes (b) acquired in SIM mode (c) XIC of PAs (500 ng/ml) in SIM mode.

LC-MS method [2]





Figure 4. LC-MS method follow BFR protocol, using a Kinetex XB-C18 column (150 x 2.1 mm, 1.7 µm). Extracted ion chromatogram of PAs (500 ng/ml) acquired in MRM mode by LCMS-8060 with ballistic gradient of 12 minutes.

Results



Method performance of LC-DMS-MS (SIM) and LC-MS (MRM)





Figure 5. Measured correction voltage (CV) at DV=800 v by infusion of diastereomeric PAs (500 ng/ml) in SIM mode. Vacuum DMS was used in scan mode (CV stepped by 0.1 volts) for optimized separation using 15% MeOH as modifier at flow rate 100 µL/min. CV values in volts for diastereomers (A) *m/z* 300 InT (-18.16), Ech (-19.46), Lyc (-18.06), InD (-16.85). (B) *m/z* 316 InT-N (-18.66), Ech-N (-22.47), Lyc-N (-21.07), InD-N (-18.36). (C) *m/z* 352 Jac (-13.34), Ret (-12.54), SeV-N (-14.44), SeC-N (-14.85). (D) *m/z* 336 SeV (-13.24) and SeC (-13.34)

QC level	Method	InT	Lyc-N	Ret	SeC
LLOQ (2 µg/kg)	LC-vDMS-MS (SIM)	Accuracy% Precision% S/N	108 1.3 205	107 7.50 17	102 2.90 17
	LC-MS (MRM)	Accuracy% Precision% S/N	102 5.40 114	102 5.40 111	102 3.50 111
	LC-vDMS-MS (SIM)	Accuracy% Precision% S/N	104 1.30 511	104 3.10 21.2	104 2.40 97.7
	LC-MS (MRM)	Accuracy% Precision% S/N	108 3.90 414	104 3.50 427	108 4.90 161.3
LQC (6 µg/kg)	LC-vDMS-MS (SIM)	Accuracy% Precision% S/N	93.4 4.80 2770	104 8.20 1374	93.4 9.00 7310
	LC-MS (MRM)	Accuracy% Precision% S/N	101 3.90 5430	104 3.50 15651	108 4.90 5317
	LC-vDMS-MS (SIM)	Accuracy% Precision% S/N	92.8 8.90 6688	111 4.20 4361	95.5 2.90 13025
	LC-MS (MRM)	Accuracy% Precision% S/N	96.5 6.10 9556	99 2.50 26153	106 2.70 8841

Figure 7. Table summarized the accuracy and precision values of PAs spiked in green tea matrix obtained at different QC levels using short LC-DMS-MS (SIM mode) and LC-MS/MS (MRM mode).

Application to quantification of PAs in tea samples




Figure 8. Quantification of selected diastereomeric pyrrolizidine alkaloids in different type of tea samples purchase in swiss supermarkets.

Conclusions

- The combined short C18 column with trap/elute LC setup hyphenated to vDMS-MS detection enables to simplify the sample preparation, higher throughput analysis (5 min) and improved selectivity compare to a classical LC-MS method.
- LC-DMS-MS enables accurate quantitation of 10 out of 14 diastereomeric PAs with higher throughput (<1h), compare to the classical LC-MS method in which 9 out of 14 diastereomeric could be separated.
- LC-vDMS-MS method in SIM mode has similar performance (linearity, accuracy and precision) compared to a reported LC-MS method in MRM.